

BD FACSCelesta™ Flow Cytometer Configuration Sheet

Blue-Violet-Yellow Green (BVYG) Laser Configuration



The BD FACSCelesta™ flow cytometer is designed to simplify the use of multicolor flow cytometry and allow researchers to benefit from new innovations in instrument and reagent technology. This platform offers multiple configurations to meet varied application needs.


The BVYG configuration of the BD FACSCelesta includes blue, violet, and yellow-green lasers, optimized to get the most out of the legacy fluorochromes and BD Horizon Brilliant™ dyes shown in the optical path diagrams. While a blue laser has long been standard in flow cytometry, the violet laser is gaining popularity as more bright violet-excited fluorochromes, such as BD Horizon Brilliant™ Violet reagents, are introduced and enthusiastically adopted. The yellow-green laser, which optimally excites phycoerythrin (PE) and bright PE-based tandem dyes as well as advanced fluorescent proteins, allows markers to be spread across more lasers, further reducing the overall compensation needed for a multicolor panel.

Multicolor flow cytometry panel design has presented challenges for researchers, such as varying marker expression, varying dye brightness, and significant emission spillover between fluorescence channels. The combination of bright, narrow-spectrum, BD Horizon Brilliant fluorochromes and sensitive optics results in panels that can readily resolve even dim populations, yet are easy to use and compensate. The system operates with BD FACSDiva™ software, a collection of convenient and easy-to-use tools for flow cytometer and application setup, data acquisition, and data analysis.

Optical path diagrams and fluorochrome support


Laser polygons show fluorochromes, mirrors, filters, and optical paths for the BD FACSCelesta BVYG configuration.

Blue laser




Other fluorochromes supported:
BD Horizon™ Fixable Viability Stain 520

Violet laser




Other fluorochromes supported:
BD Horizon™ Fixable Viability Stain 450, 510

Yellow-Green laser




Other fluorochromes supported:
BD Horizon™ Fixable Viability Stain 570, 620





8-color NK and Treg Analysis on the BD FACSCelesta BVYG Configuration.

This panel demonstrates the sensitivity and resolution of the BD FACSCelesta flow cytometer configured with the 561-nm yellow-green laser in detecting rare cell subpopulations. After normal human whole blood was washed and lysed, BD Horizon Brilliant and traditional dyes were used to identify rare NK-cell and Treg subpopulations. A. Cells were gated to select the CD3⁺ or CD3⁻ lymphocytes. B. CD3⁻ lymphocytes were analyzed to show the three subsets of NK cells based on CD16 and CD56 expression. C. Gated on the CD56^{dim}CD16⁺ NK-cell population, surface markers were used to identify CD57 and CD8 NK-cell subpopulations. D. Gated on CD3⁺ T cells, CD8⁻ T cells were identified. E. Using surface markers, Tregs were identified from the CD8⁻CD4⁺ helper T cells based on CD25 expression. F. Gated on Tregs, Treg activated and resting subsets can be resolved based on HLA-DR expression.



Expression of Fluorescent Proteins in Transfected Human Embryonic Kidney Cells.

HEK-293 cells were transfected over 24 hours with AcGFP (left) or mCherry (right), fixed with BD Cytotix™ Fixation Buffer (Cat. No. 554655), and cryopreserved for one week. After thawing, cells were washed and then analyzed by flow cytometry. Transfected cells (green or red) were compared to wild type (black) for expression of AcGFP or mCherry, respectively.

Excitation Laser	Fluorochrome	Ex _{max}	Em _{max}	Relative Brightness
Violet (405 nm)	BD Horizon Brilliant™ Violet 786 (BV786)	407 nm	786 nm	■ ■ ■ □
	BD Horizon Brilliant™ Violet 711 (BV711)	407 nm	711 nm	■ ■ ■ ■
	BD Horizon Brilliant™ Violet 650 (BV650)	407 nm	650 nm	■ ■ ■ ■
	BD Horizon Brilliant™ Violet 605 (BV605)	407 nm	602 nm	■ ■ ■ □
	BD Horizon Brilliant™ Violet 510 (BV510)	405 nm	510 nm	■ ■ ■ ■
	BD Horizon Brilliant™ Violet 421 (BV421)	407 nm	421 nm	■ ■ ■ ■
Blue (488 nm)	PerCP-Cy™5.5	482 nm	695 nm	■ ■ ■ □
	BD Horizon Brilliant™ Blue 515 (BB515)	490 nm	515 nm	■ ■ ■ ■
	FITC	494 nm	520 nm	■ ■ ■ □
Yellow-Green (561 nm)	PE-Cy™7	496 nm	785 nm	■ ■ ■ ■
	PE-Cy™5	496 nm	667 nm	■ ■ ■ ■
	BD Horizon™ PE-CF594	496 nm	612 nm	■ ■ ■ ■
	PE	496 nm	578 nm	■ ■ ■ ■

Relative Brightness Key: ■ □ □ □ Dim
 ■ ■ ■ □ Moderate
 ■ ■ ■ ■ Bright
 ■ ■ ■ ■ Brightest

Ordering Information

Description	Cat. No.
BD FACSCelesta™ Flow Cytometer, BVYG Configuration	660345
BD FACSCelesta™ High Throughput Sampler (HTS) Option	658946
BD FACSTream™ Supply System	649908
BD FACSCelesta™ Standard Workstation Bundle	660472

Class 1 Laser Product.
 For Research Use Only. Not for use in diagnostic or therapeutic procedures.
 Alexa Fluor® is a registered trademark and Pacific Blue™ is a trademark of Life Technologies Corporation.
 CF™ is a trademark of Biotium, Inc.
 Cy™ is a trademark of GE Healthcare. Cy™ dyes are subject to proprietary rights of GE Healthcare and Carnegie Mellon University, and are made and sold under license from GE Healthcare only for research and in vitro diagnostic use. Any other use requires a commercial sublicense from GE Healthcare, 800 Centennial Avenue, Piscataway, NJ 08855-1327, USA.
 BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2015 BD
 23-17601-00

