

# BD Accuri™ C6 Plus Panel Builder

| Excitation Laser Line | Fluorescence Channel  | Filter | Panel 1 | Panel 2 | Panel 3 |
|-----------------------|---|--------|---------|---------|---------|
| 488 nm                | BD Horizon™ BB515, FITC, Alexa Fluor® 488, BD Horizon™ CFSE |        |         |         |         |
|                       | PI  |        |         |         |         |
|                       | PerCP, PerCP-Cy™5.5, PE-Cy™7, 7-AAD                         |        |         |         |         |
| 640 nm                | APC, Alexa Fluor® 647                                       |        |         |         |         |

Class 1 Laser Product.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Alexa Fluor® is a registered trademark of Life Technologies Corporation.

Cy™ is a trademark of GE Healthcare. Cy™ dyes are subject to proprietary rights of GE Healthcare and Carnegie Mellon University, and are made and sold under license from GE Healthcare only for research and in vitro diagnostic use. Any other use requires a commercial sublicense from GE Healthcare, 800 Centennial Avenue, Piscataway, NJ 08855-1327, USA.

CF™ is a trademark of Biotium, Inc.

Trademarks are the property of their respective owners.

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.

23-18441-00



---

## BD Biosciences Regional Offices

### **Australia**

[bdbiosciences.com/anz](http://bdbiosciences.com/anz)

### **Canada**

[bdbiosciences.com/ca](http://bdbiosciences.com/ca)

### **China**

[bdbiosciences.com/cn](http://bdbiosciences.com/cn)

### **Europe**

[bdbiosciences.com/eu](http://bdbiosciences.com/eu)

### **India**

[bdbiosciences.com/in](http://bdbiosciences.com/in)

### **Japan**

**Nippon Becton Dickinson**  
[bd.com/jp](http://bd.com/jp)

### **Latin America/Caribbean**

[bdbiosciences.com/br](http://bdbiosciences.com/br)

### **New Zealand**

[bdbiosciences.com/anz](http://bdbiosciences.com/anz)

### **Singapore**

[bdbiosciences.com/sg](http://bdbiosciences.com/sg)

### **United States**

[bdbiosciences.com](http://bdbiosciences.com)

Office locations are available on our websites.

